



Accessories for automated slide stainers for transportation, storage, and processing of preparations on slides (applicable for cervical screening).

Vertical rack insert for 20 slides (RI-20) (Fig 1) designated for placing slides with cell samples taken from the patient, for example, during cytological screening. The capacity of the RI-20 is 20 slides. For transportation of the RI-20 in the laboratory, the special protective containers are used (see Fig. 3-4):

- plastic suitcase for 10 pcs. RI-20 with a tray for accompanying documents (capacity of at least 200 A4 sheets);
- a box made of corrugated cardboard with ten sockets (for 10 RI-20 with slides);
- a box made of corrugated cardboard with five sockets (for 5 RI-20 with slides);
- a box made of corrugated cardboard with two sockets (for 2 RI-20 with slides).

The use of the insert-racks with special containers minimizes manipulations with slides and ensures the safety of preparations.

RH-20 racks-holders (Fig. 2) are used for processing preparations in automated slide stainers FS-12-25 (FS-9-25, FS-16-COMBO, FS-16-HISTO).



Fig. 1 Rack-insert RI-20 for 20 slides



Fig. 2 Rack-holder RH-20 with RI-20

In the laboratory, the rack insert with slides (Fig. 1) are placed in rack-holders (Fig. 2) without reloading the slides, RH-20 with RI-20 are loaded into the automated slide stainer, where the preparations are stained according to a given program.

RI-20 rack Insert is supplied in protective containers (reusable or disposable). The disposable protective packaging in the form of corrugated cardboard boxes, with minor wear, can be reused.



Fig. 3 Reusable protective container-suitcase for ten rack-insert RI-20.



Fig. 4 Disposable protective packaging boxes on the 2 and 5 insert-racks RI-20.

RI-20, RH-20 are accessories of automated slide stainers FS-12-25, FS-16-COMBO, FS-16-HISTO -9-25). **Those accessories are applicable for cervical screening.**



MLT LLC: QUALITY, SIMPLICITY, AFFORDABILITY

Tehnologicheskaya street 7, 141981 Dubna, Moscow region, Russia;
+7(495) 287-81-00; www.fastainer.com; info@mlt.ru, mlt.llc@mail.ru

Version 1.2
08.11.2021

